# **EXHIBIT 8**



#### PERFORMANCE SUPPLY LLC

#### FORMAL QUOTE

TO: Ms. Ebony P. Roberson

Buyer

Office of Citywide Procurement Telephone: +1 (212) 386-0466 ERoberson@dcas.nyc.gov

1. 3M Particulate Respirator 8210, N95 160 EA/Case

Product Specification: <u>https://www.3m.com/3M/en\_US/company-us/all-3m-products/~/3M-</u>Particulate-Respirator-8210-N95-160-EA-Case/?N=5002385+3294780268&rt=rud

2. 3M Health Care Particulate Respirator 1860, N95 120 EA/Case

Product Specification: <a href="https://www.3m.com/3M/en\_US/company-us/all-3m-products/~/3M-Health-Care-Particulate-Respirator-and-Surgical-Mask-1860-N95-120-EA-Case/?N=5002385+3294795990&rt=rud">https://www.3m.com/3M/en\_US/company-us/all-3m-products/~/3M-Health-Care-Particulate-Respirator-and-Surgical-Mask-1860-N95-120-EA-Case/?N=5002385+3294795990&rt=rud</a>

3. Country of Origin: TBD

4. Headquarters: St. Paul, Minnesota, USA

5. Quantity: 8210 - 2 Million

1860 – 5 Million

6. Price: 8210 – \$6.05/mask

1860 – \$6.35/mask

- 7. Payment Terms: No upfront deposits are required. Instead, Buyer sends Purchase Order first -- and once that is accepted by the factory financial closing terms will be mutually agreed. Payment in full must be made upon receipt of Each shipment Net 15 days.
- 8. Shipping: Due to the national emergency, acceptance of the purchase order is at the full discretion of 3M and supplies are based upon availability. The N95 masks 3M can begin shipping in 2-4 weeks CIF at any of 3M plants in the USA or 3M Plants Overseas according to their manufacturing schedule. 3M chooses the plant. Order may be shipped in whole or in part.
- 9. Buyer will not terminate or reject delivery for 3M Masks after Quality Assurance inspection and products are authenticated, verified and meet 3M specifications which are indicated and agreed by the formal quote and purchase order.
- 10. Please address Purchase Order to:

Mr. Ron Romano, President Performance Supply, LLC P.O. Box 672 Englishtown, NJ 07726 Via 3BP, Inc







# **Technical Specification Sheet**

## 3M<sup>™</sup> Particulate Respirator 8210, N95

#### **Key Features**

- NIOSH approved N95 rating
- Adjustable nose clip
- Nose foam
- Ultrasonically welded headbands

#### **Material Composition**

- Straps Thermoplastic Elastomer
- Nose Clip Aluminum
- Nose foam Polyurethane
- Filter Polypropylene
- Shell Polyester
- Coverweb Polyester
- This respirator contains no components made from natural rubber latex
- Approximate weight of product: 0.35 oz.

#### **Country of Origin**

Made in the USA with globally sourced materials

#### **Use For**

- Use for solid particulates and liquid mists in concentrations not exceeding 10X PEL/OEL
- Always follow *User Instructions* and use in manners as indicated



#### Do Not Use For

- DO NOT use for gases and vapors, oil aerosols, asbestos, arsenic, cadmium, lead, 4,4-methylene dianiline (MDA), or abrasive blasting
- DO NOT use for particulate concentrations exceeding 10X PEL/OEL
- DO NOT use in any manner not indicated in the *User Instructions*

#### **Approvals and Standards**

- NIOSH approved N95 particulate respirator
- Meets NIOSH 42 CFR 84 N95
  requirements for a minimum 95%
  filtration efficiency against solid and
  liquid aerosols that do not contain oil.
- NIOSH approval number: TC-84A-0007
- Assigned Protection Factor (APF 10) per US OSHA and Canada CSA

Ordering Information								
Description	UPC	ID#	Respirators/ Box	Boxes/Case	Each/Case			
Particulate Respirator 8210	50051138464573	70-0706-1439-4	20	8	160			



#### **Time Use Limitation**

Replace the respirator when it becomes dirty, damaged, or difficult to breathe through.

#### **Shelf Life and Storage**

- 5 years from the date of manufacture
- Use By date on box in MM/YYYY format
- Store respirators in the original packaging, away from contaminated areas, dust, sunlight, extreme temperatures, excessive moisture, and damaging chemicals
- Store in temperatures between -4°F (-20°C) and +86°F (+30°C) and not exceeding 80% RH

### WARNING! A

This respirator helps reduce exposures to certain airborne contaminants. Before use, the wearer must read and understand the User Instructions provided as a part of the product packaging. Follow all local regulations. In the U.S., a written respiratory protection program must be implemented meeting all the requirements of OSHA 1910.134, including training, fit testing and medical evaluation. In Canada, CSA standard Z94.4 requirements must be met and/or requirements of the applicable jurisdiction, as appropriate. Misuse may result in sickness or death. For correct use, consult supervisor and the *User Instructions* or call 3M PSD Technical Service in USA at 1-800-243-4630 and in Canada at 1-800-267-4414.





Personal Safety Division 3M Center, Building 0235-2W-70 St. Paul, MN 55144-1000 3M.com/workersafety

#### **Acceptable Fit Test Protocols**

Fit Test P	rotocol*	Acceptable with this product?	
	Saccharin	$\boxtimes$	
Qualitative	BitrexTM	×	
Protocols	Irritant Smoke		
	Isoamyl Acetate		
Quantitative	e Protocols	$\boxtimes$	

<sup>\*</sup>Refer to OSHA 1910.134

For More Information
Technical Assistance 1-800-243-4630
Hours of Operation: M-Th 8am - 6pm, Fri 8am - 4:30 pm CST
© 3M 2018. All rights reserved.





# **Technical Specification Sheet**

## 3M<sup>™</sup> Health Care Particulate Respirator and Surgical Mask, 1860, N95

#### **Key Features**

- NIOSH approved N95 rating
- FDA cleared for use as a surgical mask
- Fluid Resistant 120 mmHg
- Flammability Rating Class I
- Adjustable nose clip
- Braided and stapled headbands

#### **Material Composition**

- Straps Braided Polyisoprene
- Staples Steel
- Nose Clip Aluminum
- Nose Foam Polyurethane Foam
- Filter Polypropylene
- Shell Polyester
- Coverweb Polypropylene
- Not made with natural rubber latex
- Approximate weight of product: 0.40 oz.

#### **Country of Origin**

Made in the USA with globally sourced materials

#### **Use For**

- Intended to be worn by operating room personnel during surgical procedures to help protect both the surgical patient and the operating room personnel from transfer of microorganisms, body fluids, and particulate material.
- Always follow *User Instructions* and use in manners as indicated



#### Do Not Use For

- DO NOT use in industrial settings
- DO NOT use for gases or vapors (i.e. anesthetic gases such as isoflurane or vapors from sterilants such as glutaraldehyde.)
- DO NOT use in any manner not indicated in the *User Instructions*

#### **Approvals and Standards**

- NIOSH approved N95 respirator
- Meets NIOSH 42 CFR 84 N95 requirements for a minimum 95% filtration efficiency against solid and liquid aerosols that do not contain oil.
- NIOSH approval number: TC-84A-0006
- FDA cleared for use as a surgical mask
- Health Canada Class I medical device
- Bacterial Filtration Efficiency F2101 >99% BFE
- Assigned Protection Factor (APF 10) per US OSHA and Canada CSA

Ordering Information								
Description	UPC	ID#	Respirators/ Box	Boxes/Case	Each/Case			
Health Care Particulate Respirator, 1860	50707387419429	70-0706-1236-4	20	6	120			



# DOCKET

# Explore Litigation Insights



Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

# **Real-Time Litigation Alerts**



Keep your litigation team up-to-date with **real-time** alerts and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

### **Advanced Docket Research**



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

## **Analytics At Your Fingertips**



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

#### API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

#### **LAW FIRMS**

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

#### **FINANCIAL INSTITUTIONS**

Litigation and bankruptcy checks for companies and debtors.

#### **E-DISCOVERY AND LEGAL VENDORS**

Sync your system to PACER to automate legal marketing.

