

'315 Patent Claims	AirFit P10 ¹
	<div data-bbox="1003 170 1465 597" data-label="Image"> </div> <p data-bbox="499 597 1087 633"><i>See Exhibit 26 (AirFit P10 Product Webpage)</i></p> <div data-bbox="787 633 1684 1019" data-label="Image"> </div> <p data-bbox="499 1026 1060 1060"><i>See Exhibit 24 (Photographs of AirFit P10).</i></p>
<p data-bbox="147 1101 472 1317">9. The mask assembly of claim 18, wherein the headgear strap comprises a flexible tube engaged with the first and second headgear extensions.</p>	<p data-bbox="499 1101 1732 1136">The AirFit P10 comprises a flexible tube engaged with the first and second headgear extensions.</p> <p data-bbox="499 1177 672 1209">For example:</p>

'315 Patent Claims



AirFit P10¹

Fitting / Mise en place / Colocación / Colocação



1. With the gray side facing out, insert each end of the headgear strap over the stabilizers through the buttonholes. Secure in place by wrapping each strap end up and around the nodes. Holding the split straps, pull the headgear from the stabilizer to ensure it is tightly secure around the nodes. Ensure the headgear strap is not twisted.

See Exhibit 25 (AirFit P10 User Guide) at 1, 2, 4.

'315 Patent Claims	AirFit P10 ¹
	 <p data-bbox="499 597 1087 634"><i>See Exhibit 26 (AirFit P10 Product Webpage)</i></p>  <p data-bbox="499 1024 1062 1060"><i>See Exhibit 24 (Photographs of AirFit P10).</i></p>