

CURRICULUM VITAE

University of Idaho

NAME: Patrick J. Hrdlicka

DATE: Sep-18

RANK: Professor

DEPARTMENT: Chemistry

OFFICE LOCATION AND CAMPUS ZIP:

University of Idaho, Department of Chemistry
875 Perimeter Drive MS 2343 (Renfrew Hall 307)
Moscow, ID 83844-2343, USA

OFFICE PHONE: 208-885-0108

FAX: 208-885-6173

EMAIL: hrdlicka@uidaho.edu

WEB: <http://www.webpages.uidaho.edu/~hrdlicka/index.htm>

DATE OF FIRST EMPLOYMENT AT UI: August 2006

DATE OF TENURE: July 2011

IMMIGRATION STATUS: Permanent resident (green card holder)

PROFESSIONAL PREPARATION:

- 2006 Ph.D (Chemistry), Univ. Southern Denmark, Odense, Denmark (supervisor: Jesper Wengel)
- 2004 M.Sc (Chemistry), Univ. Southern Denmark, Odense, Denmark (supervisor: Jesper Wengel)
- 2000 B.Sc (Chemistry), Univ. Southern Denmark, Odense, Denmark (supervisor: JIL Iversen, deceased)

EXPERIENCE:

- 2017- Professor, Dept. Chemistry, Univ. Idaho
- 2011-2017 Associate Professor, Dept. Chemistry, Univ. Idaho
- 2006-2011 Assistant Professor, Dept. Chemistry, Univ. Idaho

- 2013-2016 Affiliate Faculty, Dept. Biological Sciences, Univ. Idaho
- 2011-2014 Adjunct Faculty, Dept. Chemistry, Washington State University
- 2006-2013 Affiliate Faculty, Neuroscience Graduate Program, Univ. Idaho

- 2016- Retained as a subject matter expert on three litigations pertaining to DNA sequencing
- 2016-2017 Retained as a subject matter expert on a litigation pertaining to siRNA
- 2011-2015 Consultant for Minitube of America (now MOFA Global)
- 2011-2012 Consultant for ISIS Pharmaceuticals Inc (now Ionis Pharmaceuticals)

SCHOLARSHIP ACCOMPLISHMENTS:

Publications:

62 peer-reviewed articles and reviews (according to Google Scholar: 1350 citations, h-index = 21, i10-index = 40); eight patent applications (three issued); two book chapters; 35 published conference proceedings/abstracts; one published book review.

Peer Reviewed Journal Articles and Reviews:

[62] P. J. Hrdlicka* and S. Karmakar. 25 Years and still going strong: 2'-O-(Pyren-1-yl)methylribonucleotides - Versatile building blocks for applications in molecular biology, diagnostics and materials science. *Org. Biomol. Chem.*, **2017**, *15*, 9760-9774.

- [61] S. Karmakar, D. C. Guenther, B. C. Gibbons and P. J. Hrdlicka*. Recognition of mixed-sequence DNA using double-stranded probes with interstrand zipper arrangements of O2'-triphenylene- and coronene-functionalized RNA monomers. *Org. Biomol. Chem.*, **2017**, *15*, 9362-9371.
- [60] M. Kaura and P. J. Hrdlicka*. Efficient discrimination of single nucleotide polymorphisms (SNPs) using oligonucleotides modified with C5-pyrene-functionalized DNA and flanking Locked Nucleic Acid (LNA) monomers. *Chem. Asian J.*, **2016**, *11*, 1366-1369.
- [59] B. A. Anderson and P. J. Hrdlicka*. Merging two strategies for mixed-sequence recognition of double-stranded DNA: Pseudocomplementary Invader probes. *J. Org. Chem.*, **2016**, *81*, 3335-3346.
- [58] S. Geny, P. M. D Moreno, T. Krzywkowski, N. K. Andersen, A. J. Isse, A. M. El-Madani, C. Lou, Y. V. Pabon, O. Gissberg, B. A. Anderson, E. M. Zaghoul, R. Zain, P. J. Hrdlicka, P. T. Jørgensen, M. Nilsson, K. E. Lundin, E. B. Pedersen, J. Wengel, and C. I. E. Smith. Next generation bis-locked nucleic acids with stacking linker and 2'-glycylamino-LNA show enhanced invasion into supercoiled duplex DNA. *Nucleic Acids Res.*, **2016**, *44*, 2007-2019.
- [57] D. G. Singleton, R. Hussain, G. Siligardi, P. Kumar, P. J. Hrdlicka, N. Berova and E. Stulz*. Increased duplex stabilization in porphyrin-LNA zipper arrays with structure dependent excitonic coupling. *Org. Biomol. Chem.*, **2016**, *14*, 149-157.
- [56] D. C. Guenther, S. Karmakar and P. J. Hrdlicka*. Bulged Invader probes: Activated duplexes for mixed-sequence dsDNA recognition with improved thermodynamic and kinetic profiles. *Chem. Commun.*, **2015**, *51*, 15051-15054.
- [55] B. A. Anderson, S. Karmakar and P. J. Hrdlicka*. Mixed-sequence recognition of double-stranded DNA using enzymatically stable phosphorothioate Invader probes. *Molecules*, **2015**, *20*, 13780-13793.
- [54] B. A. Anderson and P. J. Hrdlicka*. Synthesis and characterization of oligodeoxyribonucleotides modified with 2'-thio-2'-deoxy-2'-S-(pyren-1-yl)methyluridine. *Bioorg. Med. Chem. Lett.*, **2015**, *25*, 3999-4004.
- [53] D. C. Guenther, G. H. Anderson, S. Karmakar, B. A. Anderson, B. A. Didion, W. Guo, J. P. Versteegen and P. J. Hrdlicka*. Invader probes: Harnessing the energy of intercalation to facilitate recognition of chromosomal DNA for diagnostic applications. *Chem. Sci.*, **2015**, *6*, 5006-5015.
- [52] M. Kaura and P. J. Hrdlicka*. Locked Nucleic Acid (LNA) induced effect on hybridization and fluorescence properties of oligodeoxyribonucleotides modified with nucleobase-functionalized DNA monomers. *Org. Biomol. Chem.*, **2015**, *13*, 7236-7247.
- [51] B. A. Anderson, J. J. Onley and P. J. Hrdlicka*. Recognition of double-stranded DNA using energetically activated duplexes modified with N2'-pyrene-, perylene-, or coronene-functionalized 2'-N-methyl-2'-amino-DNA monomers. *J. Org. Chem.*, **2015**, *80*, 5395-5406 (**Featured Article**).
- [50] S. Karmakar, A. S. Madsen, D. C. Guenther, B. C. Gibbons and P. J. Hrdlicka*. Recognition of double-stranded DNA using energetically activated duplexes with interstrand zippers of 1-, 2- or 4-pyrenyl-functionalized O2'-alkylated RNA monomers. *Org. Biomol. Chem.*, **2014**, *12*, 7758-7773.
- [49] B-A. F. Kengne, S. Karmakar, M. Kaura, G. Corti, V. V. R. Sai, I. B. Niraula, A. Larin, J. Hall, D. Sowell, P. J. Hrdlicka, V. Dobrokhotov and D. N. McIlroy*. Self-assembled monolayers of thiols adsorbed on Au/ZnO-functionalized silica nanosprings: Photoelectron spectroscopy-analysis and detection of vaporized explosives. *ACS Appl. Mater. Interfaces*, **2014**, *6*, 13355-13366.
- [48] D. C. Guenther, P. Kumar, B. A. Anderson and P. J. Hrdlicka*. C5-Amino acid functionalized LNA: Positively poised for antisense applications. *Chem. Commun.*, **2014**, *50*, 9007-9009.

- [47] M. Kaura, P. Kumar and P. J. Hrdlicka*. Synthesis, hybridization characteristics and fluorescence properties of oligonucleotides modified with nucleobase-functionalized LNA (Locked Nucleic Acid) adenosine and cytidine monomers. *J. Org. Chem.*, **2014**, *79*, 6256-6268.
- [46] M. Kaura, D. C. Guenther and P. J. Hrdlicka*. Carbohydrate-functionalized Locked Nucleic Acids (LNAs): Oligonucleotides with extraordinary binding affinity, target specificity and enzymatic stability. *Org. Lett.*, **2014**, *16*, 3308-3311.
- [45] P. Kumar, B. Bharal, B. A. Anderson, D. C. Guenther, M. E. Østergaard, P. K. Sharma and P. J. Hrdlicka*. C5-Alkynyl functionalized α -L-LNA: Synthesis, thermal denaturation experiments and enzymatic stability. *J. Org. Chem.*, **2014**, *79*, 5062-5073.
- [44] P. Kumar, M. E. Østergaard, B. Bharal, B. A. Anderson, D. C. Guenther, M. Kaura, D. J. Raible, P. K. Sharma and P. J. Hrdlicka*. Synthesis and biophysical properties of C5-functionalized LNA (Locked Nucleic Acid). *J. Org. Chem.*, **2014**, *79*, 5047-5061.
- [43] N. K. Andersen, B. A. Anderson, J. Wengel and P. J. Hrdlicka*. Synthesis and characterization of oligodeoxyribonucleotides modified with 2'-amino- α -L-LNA adenine monomers – High-affinity targeting of single-stranded DNA. *J. Org. Chem.*, **2013**, *78*, 12690-12702.
- [42] S. Karmakar, D. C. Guenther and P. J. Hrdlicka*. Recognition of mixed-sequence DNA duplexes: Design guidelines for Invaders based on 2'-O-(pyren-1-yl)methyl-RNA monomers. *J. Org. Chem.*, **2013**, *78*, 12040-12048.
- [41] S. P. Sau, A. S. Madsen, P. Podbevsek, N. K. Andersen, T. S. Kumar, S. Andersen, R. L. Rathje, B. A. Anderson, D. C. Guenther, S. Karmakar, P. Kumar, J. Plavec, J. Wengel and P. J. Hrdlicka*. Identification and characterization of 2nd generation Invader LNAs for mixed-sequence recognition of double-stranded DNA. *J. Org. Chem.*, **2013**, *78*, 9560-9570 (**Featured Article**).
- [40] B. Denn, S. Karmakar, D. C. Guenther and P. J. Hrdlicka*. Sandwich assay for mixed-sequence recognition of double-stranded DNA: Invader-based detection of targets specific to food pathogens. *Chem. Commun.*, **2013**, *49*, 9851-9853.
- [39] B. A. Didion, S. Karmakar, D. C. Guenther, S. Sau, J. P. Verstegen and P. J. Hrdlicka*. Invaders: Recognition of double-stranded DNA using duplexes modified with interstrand zippers of 2'-O-(pyren-1-yl)methylribonucleotides. *ChemBioChem*, **2013**, *14*, 1534-1538.
- [38] S. Karmakar and P. J. Hrdlicka*. DNA strands with alternating incorporations of LNA and 2'-O-(pyren-1-yl)methyluridine: SNP-discriminating RNA detection probes. *Chem. Sci.*, **2013**, *4*, 3447-3454.
- [37] M. Kaura, P. Kumar and P. J. Hrdlicka*. Synthesis and hybridization properties of oligonucleotides modified with 5-(1-aryl-1,2,3-triazol-4-yl)-2'-deoxyuridines. *Org. Biomol. Chem.*, **2012**, *10*, 8575-8578 (**Hot Communication**).
- [36] S. P. Sau, P. Kumar, P. K. Sharma and P. J. Hrdlicka*. Fluorescent Intercalator Displacement Replacement (FIDR) assay: Determination of relative thermodynamic and kinetic parameters in triplex formation - A case study using triplex-forming LNAs. *Nucleic Acids Res.*, **2012**, *40*, e162.
- [35] S. K. Rastogi, C. M. Gibson, J. R. Branen, D. E. Aston, A. L. Branen and P. J. Hrdlicka*. DNA detection on lateral flow test strips: Enhanced signal sensitivity using LNA-conjugated gold nanoparticles. *Chem. Commun.*, **2012**, *48*, 7714-7716.
- [34] G. Wang, M. R. Pappasani, P. Cheguru, P. J. Hrdlicka and R. A. Hill*. Gold-peptide nanoconjugate cellular uptake is modulated by serum proteins. *Nanomedicine: NBM*, **2012**, *8*, 822-832.
- [33] S. P. Sau and P. J. Hrdlicka*. C2'-Pyrene functionalized triazole linked DNA: Universal hybridization to

DNA/RNA. *J. Org. Chem.*, **2012**, *77*, 5-16 (Cover Article).

[32] S. Karmakar, B. A. Anderson, R. L. Rathje, S. Andersen, T. Jensen, P. Nielsen and P. J. Hrdlicka*. High-affinity DNA-targeting using readily accessible mimics of N2'-functionalized 2'-amino- α -L-LNA. *J. Org. Chem.*, **2011**, *76*, 7119-7131.

[31] P. Kumar, M. E. Østergaard and P. J. Hrdlicka*. Preparation of C5-functionalized Locked Nucleic Acids (LNAs). *Curr. Protocols Nucleic Acid Chem.*, **2011**, *44*, 4.43.1-4.43.22.

[30] M. E. Østergaard, P. Kumar, B. Baral, D. C. Guenther, B. A. Anderson, F. M. Ytreberg, L. Deobald, A. J. Paszczyński, P. K. Sharma and P. J. Hrdlicka*. C5-Functionalized DNA, LNA and α -L-LNA: Improved positional control of polarity-sensitive fluorophores facilitates improved SNP-typing. *Chem. Eur. J.*, **2011**, *17*, 3157-3165.

[29] M. E. Østergaard and P. J. Hrdlicka*. Pyrene-functionalized oligonucleotides and Locked Nucleic Acids (LNAs): Tools for fundamental research, diagnostics, and materials science. *Chem. Soc. Rev.*, **2011**, *40*, 5771-5788 (150+ citations as of September 18).

[28] V. V. R. Sai, D. Gangadean, I. Niraula, G. Corti, D. N. McIlroy, D. E. Aston, J. R. Branen and P. J. Hrdlicka*. Silica nanosprings coated with noble metal nanoparticles: Highly active SERS substrates. *J. Phys. Chem. C.*, **2011**, *115*, 453-459.

[27] M. E. Østergaard, J. Maiti, B. R. Babu, J. Wengel and P. J. Hrdlicka*. Glowing LNA as nucleic acid detection probes – Nucleobase analogs and insights into probe design. *Bioorg. Med. Chem. Lett.*, **2010**, *20*, 7265-7268.

[26] R. N. Veedu*, H. R. Burri, P. Kumar, P. K. Sharma, P. J. Hrdlicka, B. Vester and J. Wengel. Polymerase-directed synthesis of C5-ethynyl Locked Nucleic Acids. *Bioorg. Med. Chem. Lett.*, **2010**, *22*, 6565-6568.

[25] M. E. Østergaard, P. Cheguru, M. R. Papasani, R. A. Hill* and P. J. Hrdlicka*. Glowing LNA – Brightly fluorescent probes for detection of nucleic acids in cells. *J. Am. Chem. Soc.* **2010**, *132*, 14221-14228.

[24] M. E. Østergaard, D. C. Guenther, P. Kumar, B. Baral, L. Deobald, A. J. Paszczyński, P. K. Sharma and P. J. Hrdlicka*. Pyrene-functionalized triazole-linked 2'-deoxyuridines - Probes for discrimination of single nucleotide polymorphisms (SNPs). *Chem. Commun.*, **2010**, 4929-4931.

[23] S. P. Sau, T. S. Kumar and P. J. Hrdlicka*. Invader LNA – Efficient targeting of short DNA duplexes. *Org. Biomol. Chem.*, **2010**, *8*, 2028-2036.

[22] M. E. Østergaard, P. Kumar, B. Baral, D. J. Raible, T. S. Kumar, B. A. Anderson, D. C. Guenther, L. Deobald, A. J. Paszczyński, P. K. Sharma and P. J. Hrdlicka*. C5-Functionalized LNA: Unparalleled hybridization properties and enzymatic stability. *ChemBioChem*, **2009**, *10*, 2740-2743.

[21] S. P. Sau, P. Kumar, B. A. Anderson, M. E. Østergaard, L. Deobald, A. J. Paszczyński, P. K. Sharma and P. J. Hrdlicka*. Optimized DNA-targeting using triplex forming C5-alkynyl functionalized LNA. *Chem. Commun.*, **2009**, 6756-6758.

[20] S. K. Rastogi, N. N. Mishra, M. E. Østergaard, E. Cameron, B. Finaloski, P. J. Hrdlicka and W. C. Maki*. Hybridization study of PNA-DNA in liquid-solid phase for biosensor application. *Anal. Lett.*, **2009**, *42*, 2485-2495.

[19] B. Baral, P. Kumar, B. A. Anderson, M. E. Østergaard, P. K. Sharma and P. J. Hrdlicka*. Optimized synthesis of [3-¹⁵N]-labeled uridine phosphoramidites. *Tetrahedron Lett.*, **2009**, *50*, 5850-5852.

[18] T. S. Kumar, M. E. Østergaard, P. K. Sharma, P. Nielsen, J. Wengel and P. J. Hrdlicka*. Parallel RNA-strand recognition by 2'-amino- β -L-LNA. *Bioorg. Med. Chem. Lett.*, **2009**, *19*, 2396-2399.

[17] J. W. Rigoli, M. E. Østergaard, K. M. Canady, D. C. Guenther and P. J. Hrdlicka*. Selective deacylation

- [16] T. S. Kumar, A. S. Madsen, M. E. Østergaard, S. P. Sau, J. Wengel and P. J. Hrdlicka*. Functionalized 2'-amino- α -L-LNA - Directed positioning of intercalators for DNA targeting. *J. Org. Chem.* **2009**, *74*, 1070-1081.
- [15] T. S. Kumar, P. Kumar, P. K. Sharma and P. J. Hrdlicka*. Optimized synthesis of LNA uracil nucleosides. *Tetrahedron Lett.*, **2008**, *49*, 7168-7170.
- [14] T. S. Kumar, A. S. Madsen, M. E. Østergaard, J. Wengel* and P. J. Hrdlicka*. Nucleic acid structural engineering using pyrene functionalized 2'-amino- α -L-LNA monomers and abasic sites. *J. Org. Chem.* **2008**, *73*, 7060-7066.
- [13] B. Vester, A. M. Boel, S. Lobedan, B. R. Babu, M. Raunkjær, D. Lindegaard, Raunak, P. J. Hrdlicka, T. Højland, P. K. Sharma, T. S. Kumar, P. Nielsen and J. Wengel*. Chemically modified oligonucleotides with retained RNaseH response. *Bioorg. Med. Chem. Lett.* **2008**, *18*, 2296-2300.
- [12] T. Umemoto, P. J. Hrdlicka, B. R. Babu and J. Wengel*. Sensitive SNP dual-probe assays based on pyrene-functionalized 2'-amino-LNA: Lessons to be learned. *ChemBioChem*, **2007**, *8*, 2240-2248.
- [11] T. S. Kumar, J. Wengel and P. J. Hrdlicka*. 2'-N-(pyren-1-yl)acetyl-2'-amino- α -L-LNA: Synthesis and detection of single nucleotide mismatches in DNA and RNA targets. *ChemBioChem*, **2007**, *8*, 1122-1125.
- [10] A. S. Madsen, P. J. Hrdlicka, T. S. Kumar and J. Wengel*. Synthesis, nucleic acid hybridization properties and molecular modeling studies of conformationally restricted 3'-O,4'-C-methylene-linked α -L-ribonucleotides. *Carbohydr. Res.* **2006**, *341*, 1398-1407.
- [9] T. S. Kumar, A. S. Madsen, J. Wengel and P. J. Hrdlicka*. Synthesis and hybridization studies of 2'-amino- α -L-LNA and a tetracyclic 'Locked LNA'. *J. Org. Chem.* **2006**, *71*, 4188-4201.
- [8] P. J. Hrdlicka, T. S. Kumar and J. Wengel*. Synthesis and thermal denaturation studies of conformationally restricted 3'-C-ethynyl-3'-O,4'-C-methylene-ribonucleotides. *Eur. J. Org. Chem.* **2005**, 5184-5188.
- [7] P. J. Hrdlicka, B. R. Babu, M. D. Sørensen, N. Harrit and J. Wengel*. Multilabeled pyrene-functionalized 2'-amino-LNA probes for nucleic acid detection in homogenous fluorescence assays. *J. Am. Chem. Soc.* **2005**, *127*, 13293-13299 (**125 citations as of April 2018**).
- [6] P. J. Hrdlicka, T. S. Kumar and J. Wengel*. Targeting of mixed sequence double-stranded DNA using pyrene-functionalized 2'-amino- α -L-LNA. *Chem. Commun.* **2005**, 4279-4281.
- [5] B. R. Babu, P. J. Hrdlicka, C. J. McKenzie and J. Wengel*. Optimized DNA targeting using *N,N*-bis(2-pyridylmethyl)- β -alanyl 2'-Amino-LNA. *Chem. Commun.* **2005**, 1705-1707.
- [4] P. J. Hrdlicka, N. K. Andersen, J. S. Jepsen, F. G. Hansen, K. F. Haselmann, C. Nielsen and J. Wengel*. Synthesis and biological evaluation of branched and conformationally restricted analogs of the anticancer compounds 3'-C-ethynyluridine (EUrd) and 3'-C-ethynylcytidine (ECyd). *Bioorg. Med. Chem.* **2005**, *13*, 2597-2621.
- [3] P. J. Hrdlicka, J. S. Jepsen, C. Nielsen and J. Wengel*. Synthesis and biological evaluation of nucleobase-modified analogs of the anticancer compounds 3'-C-ethynyluridine (EUrd) and 3'-C-ethynylcytidine (ECyd). *Bioorg. Med. Chem.* **2005**, *13*, 1249-1260.
- [2] P. J. Hrdlicka, B. R. Babu, M. D. Sørensen and J. Wengel*. Interstrand communication between 2'-N-(pyren-1-yl)methyl-2'-amino LNA monomers in nucleic acid duplexes: Directional control and signaling of full complementarity. *Chem. Commun.* **2004**, 1478-1479.

Explore Litigation Insights

Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

Real-Time Litigation Alerts



Keep your litigation team up-to-date with **real-time alerts** and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

Advanced Docket Research



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

Analytics At Your Fingertips



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

LAW FIRMS

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

FINANCIAL INSTITUTIONS

Litigation and bankruptcy checks for companies and debtors.

E-DISCOVERY AND LEGAL VENDORS

Sync your system to PACER to automate legal marketing.